



# Repaira<sup>®</sup>

Polynucleotide Hydrogel,  
Skin's Reset Button !

An Innovation by **Global Experts**

Beyond Skin Treatments – A New Era of Intracellular & Extracellular Regeneration. Support, Repair, and Restructuring Skin for Enhanced Texture and Appearance.

**Repaira®** is the next-generation **PDRN** Exosomal DNA, designed to repair, rebuild, and restore skin at both intracellular and extracellular levels.

Powered by cutting-edge bio-engineered regenerative Growth Factors, unique technologies, and potent cellular stimulants, **Repaira®** delivers unmatched skin recovery and restructuring.



## Repaira® Unique Technologies

**GDN Exo™ in Repaira:** The Evolution of Exosomal DNA Repair Technology.

Beyond Traditional Exosomes – A Bio-Engineered Growth DNA Complex, a fusion of Growth-DNA Nano-Complex (GDN) and Exosomal therapy for precision cellular restoration.

### GDN Exo™

- Synthetic Exosome
- **PDRN (Salmon DNA)**
- Bio-engineered GF

### CAAC™:

Collagen Amino Acid Complex


**Hydro-Lattice Encapsulation™:**  
Self-Crosslinked HA

**Carnosine**  
Crucial for Repairing Skin Barrier  
And protecting Collagen & Elastin

**Recovery Plus**  
Tranexamic Acid  
Citric Acid  
Centella Asiatica  
Arginine  
Acetyl Hexapeptide-8



**Hydro-Lattice Encapsulation™:**  
Long-lasting release and precision delivery of actives.



**Self-Crosslinked HA:** Delivers superior hydration, elasticity, and firmness.



**CAAC™ Technology:** Promotes collagen renewal and structural repair.

**Repara® is clinically designed for:** (Face - Neck - Under Eyes),  
To Repair and Restructure Skin at the Intracellular & Extracellular Level.

## How Repara® Works ? Mechanism of Action

### Intracellular Resetting:

- DNA Repair.
- Cellular Activation.
- Mitochondrial Optimization.
- Neuromuscular Relaxation.

### Extracellular Matrix Remodeling:

- Collagen Biosynthesis.
- Skin Barrier Repair.
- Antioxidant Protection.
- Pigmentation & Vascular Control.

### Repara® Occasional Rejuvenation:


- 1 session 10-14 days before occasion.

### Repara® Treatment Protocol

- Intensive Reparation Phase: 2-4 Repara® sessions, spaced 2 weeks apart.
- Collagen Maturation Phase: 1 Gloform® Generate session after 1 month.
- Maintenance & Prevention: 1 Repara® / Gloform® session every 6 months for sustained ECM stability.



## Grid Mapping Method





# Repaira® Skin's Reset Button!

**Repaira®** is not just another PDRN product, it's a revolution in skin regeneration.

## Why Repaira®?

### The Science Behind the Name

The name **Repaira®** is derived from "Repair", representing its superior ability to restore, heal, and rejuvenate skin beyond conventional injectables. Unlike traditional skin boosters, **Repaira®** is not just about regeneration—it is a true skin RESET mechanism.

1. Resets Cellular Communication → Acetyl Hexapeptide-8 rewires neuromuscular micro-signaling, improving skin tone and reducing dynamic wrinkles.
2. Resets ECM Integrity → High-dose PDRN within exosomes accelerates collagen cross-linking & tissue stabilization, leading to long-term regenerative benefits.
3. Resets Aging Damage → Carnosine prevents glycation and oxidative degradation, ensuring collagen fibers remain healthy, elastic, and functional over time.
4. Resets Hydration & Barrier Function → Hydro-Lattice Encapsulation™ ensures sustained hydration and protection, reinforcing epidermal resilience.



## Science Meets Responsibility

- We prioritize biodegradable materials, sustainable sourcing, and ethical production methods, ensuring minimal environmental impact.
- Our cruelty-free and responsibly sourced ingredients reflect our dedication to ethical beauty and patient well-being.

### CONTACT US

#### GLOBAL EXPERTS SRL

[www.globalexpertsgroup.com](http://www.globalexpertsgroup.com)

[info@globalexpertsgroup.com](mailto:info@globalexpertsgroup.com)

Phone: +39-3447727704

